



HBV Real-TM Quant

new version

Handbook

Real Time Kit for the Quantitative detection
of Hepatitis B Virus in human plasma

REF V5-100/2 FRT

REF TV5-100/2 FRT

REF TV5-100/2 FRT C

 **100**

NAME

HBV Real-TM Quant

INTRODUCTION

Hepatitis B virus (HBV) is a member of the Hepadnavirus family. The virus particle, (virion) consists of an outer lipid envelope and an icosahedral nucleocapsid core composed of protein. These virions are 42 nm in diameter and are sometimes referred to as "Dane particles". The nucleocapsid encloses the viral DNA and a DNA polymerase that has reverse transcriptase activity. The outer envelope contains embedded proteins that are involved in viral binding of, and entry into, susceptible cells. The virus is one of the smallest enveloped animal viruses, but pleomorphic forms exist, including filamentous and spherical bodies lacking a core. These particles are not infectious and are composed of the lipid and protein that forms part of the surface of the virion, which is called the surface antigen (HBsAg), and is produced in excess during the life cycle of the virus.

Acute infection with hepatitis B virus is associated with acute viral hepatitis – an illness that begins with general ill-health, loss of appetite, nausea, vomiting, body aches, mild fever, and dark urine, and then progresses to development of jaundice. It has been noted that itchy skin has been an indication as a possible symptom of all hepatitis virus types. The illness lasts for a few weeks and then gradually improves in most affected people. A few people may have more severe liver disease (fulminant hepatic failure), and may die as a result. The infection may be entirely asymptomatic and may go unrecognized.^[16]

Chronic infection with hepatitis B virus either may be asymptomatic or may be associated with a chronic inflammation of the liver (chronic hepatitis), leading to cirrhosis over a period of several years. This type of infection dramatically increases the incidence of hepatocellular carcinoma (liver cancer). Chronic carriers are encouraged to avoid consuming alcohol as it increases their risk for cirrhosis and liver cancer. Approximately 300 million individuals are chronically infected with hepatitis B virus in the world. Enzyme-linked immunosorbent assay (ELISA) is still a main detection method for HBV infection, but ELISA result can neither efficiently reflect serum viral load or hepatitis activity nor monitor the efficacy of antiviral treatments. Currently, polymerase chain reaction (PCR) assay has been widely used for monitoring HBV load. HBV DNA monitoring has become an important tool to identify individuals with high viral replication, to monitor patients on therapy, and to predict whether antiviral therapy is successful. For example, with the introduction of new antiviral agents like lamivudine, close monitoring of patients has become increasingly important due to the occurrence of antiviral drug-resistant virus strains or the presence of flares after withdrawal of antiviral therapy.

INTENDED USE

HBV Real-TM Quant is a Real-Time test for the Quantitative detection of Hepatitis B Virus in human plasma and simultaneous detection of a HBV-specific Internal Control (IC), by dual color detection.

PRINCIPLE OF ASSAY

HBV Real-TM Quant is a Real-Time test for the Quantitative detection of Hepatitis B Virus in human plasma. HBV DNA is extracted from plasma, amplified using Real Time Amplification and detected using fluorescent reporter dye probes specific for HBV or HBV IC. Internal Control (IC) serves as an amplification control for each individually processed specimen and to identify possible inhibition. IC is

detected in a channel other than the HBV DNA. Monitoring the fluorescence intensities during Real Time allows the detection of the accumulating product without having to reopen the reaction tube after the amplification.

MATERIALS PROVIDED

REF V5-100/2 FRT

| HBV Real-TM Quant – PCR kit | Volume, ml | Quantity |
|---|-------------------|-----------------|
| PCR-mix-1-TM HBV | 0.3 | 4 tubes |
| PCR-mix-2-FRT | 0.2 | 4 tubes |
| Polymerase (TaqF) | 0.02 | 4 tubes |
| QS1 HBV¹ | 0.1 | 4 tubes |
| QS2 HBV¹ | 0.1 | 4 tubes |
| Negative Control (C-)² | 1,2 ml | 4 tubes |
| TE-buffer³ | 1.2 | 4 tubes |
| Positive Control-1-HBV⁴ | 0.06 | 4 tubes |
| Positive Control-2-HBV⁴ | 0.06 | 4 tubes |
| Internal Control STI-87 (IC)⁵ | 0.28 | 4 tubes |

1 must be used as Calibration standards. Prepare two (2) tubes for each standard during the amplification procedure. Concentrations are specific for every lot (see DataCard).

2 must be used in the DNA extraction procedure as Negative Control of Extraction.

3 must be used in the amplification procedure as Negative Control of Amplification.

4 must be used in the DNA extraction procedure as Positive Controls of Extraction. If the extraction kit starts from 100 µl, add 90 µl of Negative Control (C-) and 10 µl of Positive Control-1-HBV and Positive Control-2-HBV to the corresponding tubes labeled Cpos1 and Cpos2. For example, if it starts from 400 µl add 390 µl of Negative Control (C-) and 10 µl of Positive Control-1-HBV or Positive Control-2-HBV respectively. Controls concentration range is specific for every lot (see DataCard);

5 add 10 µl of Internal Control during the DNA extraction procedure directly to the sample/lysis mixture.

| HBV Real-TM Quant – PCR kit | Volume, ml | Quantity |
|---|-------------------|-----------------|
| PCR-mix-1-TM HBV | 0.3 | 4 tubes |
| PCR-mix-2-FRT | 0.2 | 4 tubes |
| Polymerase (TaqF) | 0.02 | 4 tubes |
| QS1 HBV ¹ | 0.1 | 4 tubes |
| QS2 HBV ¹ | 0.1 | 4 tubes |
| Negative Control (C-) ² | 1,2 ml | 4 tubes |
| TE-buffer ³ | 1.2 | 4 tubes |
| Positive Control-1-HBV ⁴ | 0.06 | 4 tubes |
| Positive Control-2-HBV ⁴ | 0.06 | 4 tubes |
| Internal Control STI-87 (IC) ⁵ | 0.28 | 4 tubes |
| Ribo-Sorb – Extraction kit | Volume, ml | Quantity |
| Lysis Solution | 22.5 | 2 tubes |
| Washing Solution | 20 | 2 tubes |
| Sorbent | 1.25 | 2 tubes |
| RNA-eluent | 0.5 | 10 tubes |

- 1 must be used as Calibration standards. Prepare two (2) tubes for each standard during the amplification procedure. Concentrations are specific for every lot (see DataCard).
- 2 must be used in the DNA extraction procedure as Negative Control of Extraction.
- 3 must be used in the amplification procedure as Negative Control of Amplification.
- 4 must be used in the DNA extraction procedure as Positive Controls of Extraction. If the extraction kit starts from 100 µl, add 90 µl of Negative Control (C-) and 10 µl of Positive Control-1-HBV and Positive Control-2-HBV to the corresponding tubes labeled Cpos1 and Cpos2. For example, if it starts from 400 µl add 390 µl of Negative Control (C-) and 10 µl of Positive Control-1-HBV or Positive Control-2-HBV respectively. Controls concentration range is specific for every lot (see DataCard);
- 5 add 10 µl of Internal Control during the DNA extraction procedure directly to the sample/lysis mixture.

| HBV Real-TM Quant – PCR kit | Volume, ml | Quantity |
|--|-------------------|-----------------|
| PCR-mix-1-TM HBV | 0.3 | 4 tubes |
| PCR-mix-2-FRT | 0.2 | 4 tubes |
| Polymerase (TaqF) | 0.02 | 4 tubes |
| QS1 HBV ¹ | 0.1 | 4 tubes |
| QS2 HBV ¹ | 0.1 | 4 tubes |
| Negative Control (C-) ² | 1,2 ml | 4 tubes |
| TE-buffer ³ | 1.2 | 4 tubes |
| Positive Control-1-HBV ⁴ | 0.06 | 4 tubes |
| Positive Control-2-HBV ⁴ | 0.06 | 4 tubes |
| Internal Control STI-87 (IC) ⁵ | 0.28 | 4 tubes |
| Ribo-Virus – Extraction kit | Volume, ml | Quantity |
| Buffer RAV1 | 35 | 2 tubes |
| Buffer RAW | 30 | 2 tubes |
| Buffer RAV3 (concentrate) | 12 | 2 tubes |
| Buffer RE | 13 | 2 tubes |
| Rnase-free H ₂ O | 13 | 2 tubes |
| Carrier RNA (lyophilized) | 1 (mg) | 2 tubes |
| Proteinase K | 50 (mg) | 2 tubes |
| Proteinase buffer | 8 | 1 tube |
| Ribo Virus columns with collecting tubes (2ml) | / | 100 |
| Collecting tubes (2ml) | / | 8 x 50 |

- 1 must be used as Calibration standards. Prepare two (2) tubes for each standard during the amplification procedure. Concentrations are specific for every lot (see DataCard).
- 2 must be used in the DNA extraction procedure as Negative Control of Extraction.
- 3 must be used in the amplification procedure as Negative Control of Amplification.
- 4 must be used in the DNA extraction procedure as Positive Controls of Extraction. If the extraction kit starts from 100 µl, add 90 µl of Negative Control (C-) and 10 µl of Positive Control-1-HBV and Positive Control-2-HBV to the corresponding tubes labeled Cpos1 and Cpos2. For example, if it starts from 400 µl add 390 µl of Negative Control (C-) and 10 µl of Positive Control-1-HBV or Positive Control-2-HBV respectively. Controls concentration range is specific for every lot (see DataCard);
- 5 add 10 µl of Internal Control during the DNA extraction procedure directly to the sample/lysis mixture.


MATERIALS REQUIRED BUT NOT PROVIDED

- DNA isolation kit (only ref. V5-100/2 FRT)
- Biological cabinet
- Desktop microcentrifuge for “Eppendorf” type tubes (RCF max. 16,000 x g); Eppendorf 5415D or equivalent
- 60°C ± 2°C dry heat block
- Vortex mixer
- Pipettors (capacity 5-40 µl; 40-200 µl; 200-1000 µl) with aerosol barrier
- 1,5 ml polypropylene sterile tubes (Sarstedt, QSP, Eppendorf)
- Disposable gloves, powderless
- Tube racks
- Biohazard waste container
- 70% Ethanol (freshly prepared mixture of reagent grade 96% ethanol and distilled water)
- Acetone
- Refrigerator
- Real Time Thermal cycler
- Freezer

WARNINGS AND PRECAUTIONS



For Research Use Only

1.  Lysis Solution contains guanidine thiocyanate. Guanidine thiocyanate is harmful if inhaled, or comes in contact with skin or if swallowed. Contact with acid releases toxic gas. (Xn; R: 20/21/22-36/37/38; S: 36/37/39) - (only ref. TV5-100/2FRT – TV5-100/2FRT C)
2. Wear disposable gloves, laboratory coats and eye protection when handling specimens and reagents. Thoroughly wash hands afterward.
3. Use routine laboratory precautions. Do not eat, drink, smoke, apply cosmetics, or handle contact lenses in laboratory work areas. Do not pipette by mouth.
4. Do not use a kit after its expiration date.
5. Do not mix reagents from different kits.
6. Dispose all specimens and unused reagents in accordance with local regulations.
7. Heparin has been shown to inhibit reaction. The use of heparinized specimens is not recommended.
8. Avoid repeated thawing and freezing of the reagents, this may reduce the sensitivity of the test.
9. Once the reagents have been thawed, vortex and centrifuge briefly the tubes.
10. Prepare quickly the Reaction mix.
11. Specimens may be infectious. Use Universal Precautions when performing the assay.
12. Specimens and controls should be prepared in a laminar flow hood.
13. Handle all materials containing specimens or controls according to Good Laboratory Practices in order to prevent cross-contamination of specimens or controls.
14. Clean and disinfect all spills of specimens or reagents using a disinfectant such as 0,5% sodium hypochlorite, or other suitable disinfectant. Follow by wiping down the surface with 70% ethanol.
15. Avoid contact of specimens and reagents with the skin, eyes and mucous membranes. If these solutions come into contact, rinse immediately with water and seek medical advice immediately.
16. Material Safety Data Sheets (MSDS) are available on request.
17. Use of this product should be limited to personnel trained in the techniques of amplification.
18. Workflow in the laboratory must proceed in a uni-directional manner, beginning in the Extraction Area and moving to the Amplification Area. Do not return samples, equipment and reagents in the area where you performed previous step. Personnel should be using proper anti-contamination safeguards when moving between areas.

STORAGE INSTRUCTIONS

HBV Real-TM Quant must be stored at temperature $\leq -20^{\circ}\text{C}$. The kit can be shipped at $2-8^{\circ}\text{C}$ for 3-4 days but should be stored at $\leq -20^{\circ}\text{C}$ immediately on receipt.

STABILITY

HBV Real-TM Quant Test is stable up to the expiration date indicated on the kit label. The product will maintain performance through the control date printed on the label. Exposure to light, heat or humidity may affect the shelf life of some of the kit components and should be avoided. Repeated thawing and freezing of these reagents should be avoided, as this may reduce the sensitivity. Components stored

under conditions other than those stated on the labels may not perform properly and may adversely affect the assay results.

SAMPLE COLLECTION, STORAGE AND TRANSPORT

Note: Handle all specimens as if they are potentially infectious agents.

1. EDTA tubes may be used with the **HBV Real-TM Quant**. Follow sample tube manufacturer's instructions.
2. Whole blood collected in EDTA should be separated into plasma and cellular components by centrifugation at 800-1600 x g for 20 min within six hours. The isolated plasma has to be transferred into a sterile polypropylene tube. Plasma may be stored at 2-8°C for an additional 3 days. Alternatively, plasma may be stored at -18°C for up to one month or 1 year when stored at -70°C.
3. Do not freeze whole blood.
4. Specimens anti-coagulated with heparin are unsuitable for this test.
5. Thaw frozen specimens at room temperature before using.
6. Whole blood must be transported at 2-25°C and processed within 6 hours of collection. Plasma may be transported at 2-8°C or frozen.
7. Transportation of clinical specimens must comply with country, federal, state and local regulations for the transport of etiologic agents.

DNA ISOLATION

The following isolation kits are recommended:

- ⇒ **DNA/RNA-Prep** (Sacace, [REF K-2-9](#))
- ⇒ **Ribo-Sorb** (Sacace, [REF K-2-1](#))
- ⇒ **Magno Virus** (Sacace, [REF K-2-16/1000](#))
- ⇒ **Ribo Virus**– spin column extraction kit (Sacace, [REF K-2-C](#))
- ⇒ **SaMag Viral Nucleic Acids Extraction kit** (Sacace, [REF SM003](#))

Please carry out the DNA extraction according to the manufacturer's instructions. Add 10 µl of Internal Control (IC) during the RNA isolation procedure directly to the sample/lysis mixture.

PCR PREPARATION

The total reaction volume is 25 µl, the volume of DNA sample is 12.5 µl.

1. Thaw one set of reagents, vortex and centrifuge briefly the tubes.
2. Prepare requested quantity of reaction tubes including 3 extraction controls, negative amplification control and 4 standards;
3. Prepare in a new tube for N samples (including controls and standards) the **Reaction Mix** as follow:
 - **7.5 µl *N** of **PCR-mix-1-TM HBV**,
 - **5 µl *N** of **PCR-mix-2-FRT**,
 - **0.5 µl *N** of **Polymerase (TaqF)**,

Vortex thoroughly and centrifuge briefly.

4. Add **12.5 µl** of **Reaction Mix** into each tube.
5. Add **12.5 µl** of **extracted DNA** sample to the appropriate tubes with Reaction Mix.
6. Prepare for each run 4 standards:

- add **12.5 µl** of **QS1 HBV** into two (2) tubes each labelled QS1 HBV. Mix by pipetting;
 - add **12.5 µl** of **QS2 HBV** into two (2) tubes each labelled QS2 HBV. Mix by pipetting;
7. To rule out any possible contamination, run an additional negative control of amplification reaction by adding **12.5 µl** of **TE-buffer** to the appropriate tube with Reaction Mix;
- Insert the tubes in the thermalcycler and program the instrument.

Amplification program

Create a temperature profile on your instrument as follows:

| Step | Rotor-type Instruments ¹ | | | Plate- or modular type Instruments ² | | |
|------|-------------------------------------|---|---------|---|---|---------|
| | Temperature, °C | Time | Repeats | Temperature, °C | Time | Repeats |
| 1 | 95 | 15 min | 1 | 95 | 15 min | 1 |
| 2 | 95 | 5 s | 5 | 95 | 5 s | 5 |
| | 60 | 20 s | | 60 | 20 s | |
| | 72 | 15 s | | 72 | 15 s | |
| 3 | 95 | 5 s | 40 | 95 | 5 s | 40 |
| | 60 | 20 s <i>fluorescent signal detection</i> | | 60 | 30 s <i>fluorescent signal detection</i> | |
| | 72 | 15 s | | 72 | 15 s | |

¹ For example Rotor-Gene™ 3000/6000/Q (Corbett Research, Qiagen)

² For example, SaCycler-96™ (Sacace), CFX-96/iQ5™ (BioRad); Mx3000/3005P™ (Agilent), ABI® 7300/7500/StepOne Real Time PCR (Applied Biosystems), SmartCycler® (Cepheid)

Fluorescence is detected at the 3th step of Cycling 2 stage (60 °C) in FAM/Green and JOE/Yellow/Hex/Cy3 fluorescence channels.

Internal control is detected on the FAM (Green) channel, HBV DNA is detected on the JOE(Yellow)/HEX/Cy3 channel.

INSTRUMENT SETTINGS

RotorGene (Corbett Research, Qiagen):

| Channel | Calibrate/Gain Optimisation... | Threshold | More Settings/ Outlier Removal | Slope Correct |
|------------|--------------------------------|-----------|--------------------------------|---------------|
| FAM/Green | from 5 Fl to 10 Fl | 0.03 | 10 % | On |
| JOE/Yellow | from 5 Fl to 10 Fl | 0.03 | 10 % | On |

Plate-type instruments like SaCycler-96* (Sacace), CFX-96 / iQ5 (BioRad), ABI 7300/7500/StepOne (Applied Biosystem):

| Channel | Threshold |
|---------|--|
| FAM | For each channel at a time set the threshold line at the level of 10-20 % of maximum fluorescence obtained for the Calibrator (QS1-HBV) in the last amplification cycle. |
| HEX | |

* for SaCycler-96 instrument set "Criterion of the PCR positive result" to 70%, select "Multiplex Detection" as analysis type and "Threshold (Ct)" as Method.

SmartCycler (Cepheid):

| Channel | Threshold |
|---------|-----------|
| FAM | 20 |
| Cy3 | 30 |

PROTOCOL AND DATA ANALYSIS:

Principle of interpretation is the following:

- The signal of the Internal Control DNA amplification product is detected in the FAM channel.
- The signal of the HBV DNA amplification product is detected in the JOE channel.

The results are interpreted by the presence (or absence) of the intercept between the fluorescence curve and the threshold line which determines the presence (or absence) of the Ct values for the sample.

Based on the Ct values and on the specified values of the calibrators, QS1 HBV and QS2 HBV, the calibration line will give the values for the number of HBV DNA copies (JOE channel) and for the number of Internal Control DNA copies (FAM channel) in a PCR sample.

The retrieved values are used for the HBV DNA concentration calculation, using the following formula:

$$\frac{\text{HBV DNA copies per PCR-sample}}{\text{IC DNA copies per PCR-sample}} \times \text{coefficient A} \times \text{coefficient B} = \text{IU/ml of plasma}$$

$$\text{Coefficient A} = \frac{100}{\text{extraction volume, } \mu\text{l}}$$



Coefficient A = 1 when calculating Positive Control-1-HBV and Positive Control-2-HBV concentrations.

For example, if using MagnoVirus extraction kit which starts from 1000 μl of plasma, coefficient A will be $100/1000 = 0,1$

Coefficient B is specified in the DataCard provided with the PCR kit and is specific for each lot. It cannot be used with PCR kits of different lots.



If the result is greater than 100,000,000 IU/ml, then it is interpreted as the greater than 100,000,000 IU/ml. If the obtained value is higher than the linear range, then the sample may be re-tested after 10x dilution; the produced result must be multiplied by 10.

If the result is less than 30 IU/ml (extraction from 1 ml), then it is interpreted as the less than 30 IU/ml.

The result of the analysis is considered reliable only if the results obtained for Positive and Negative Controls of amplification as well as for the Positive and Negative Control of extraction are correct (see table: Results for controls).

Results for controls and calibrators

| Control | Stage for control | FAM/Green (IC) | JOE/Yellow/HEX (HBV) |
|--------------------------|-------------------|----------------|----------------------|
| C- | DNA extraction | + | - |
| PCE (C+) | DNA extraction | + | + |
| NCA | Amplification | - | - |
| QS1 HBV / QS2 HBV | Amplification | + | + |



Positive Control-1-HBV and Positive Control-2-HBV results must fall in the expected range of quantitation indicated in the DataCard.



Calibration curve is valid if R² value is higher than 0,98.

QUALITY CONTROL PROCEDURE

A defined quantity of Internal Control (IC) is introduced into each sample and control at the beginning of sample preparation procedure in order to control the extraction process of each individual sample and to identify possible reaction inhibition.

A negative control of extraction (NCE), negative amplification control (NCA), positive control of extraction (C+), positive amplification control (QS2 HBV) are required for every run to verify that the specimen preparation, the amplification and the detection steps are performed correctly.

If the controls are out of their expected range (see table Results for Controls), all of the specimens and controls from that run must be processed beginning from the sample preparation step.

PERFORMANCE CHARACTERISTICS

Analytical specificity

The analytical specificity of **HBV Real-TM Quant** PCR kit is ensured by selection of specific primers and probes and stringent reaction conditions. The primers and probes were tested for possible homologies to all sequences deposited in gene banks by sequence comparison analysis as well as with genomic DNA/RNA of the following organisms and viruses: hepatitis A virus; hepatitis C virus; hepatitis D virus; human immunodeficiency virus; cytomegalovirus; Epstein-Barr virus; herpes simplex virus types 1 and 2; *varicella-zoster* virus; human herpes virus types 6 and 8; parvovirus B19; tick-borne encephalitis virus; West Nile encephalitis; adenovirus types 2, 3, and 7; *Escherichia coli*; *Staphylococcus aureus*; *Streptococcus pyogenes*, *S.agalactiae*; and *Homo sapiens*. Cross-reactions for the above-mentioned organisms and viruses have not been detected.

Analytical sensitivity

The kit **HBV Real-TM Quant** allows to detect *HBV* DNA in 100% of the tests with a sensitivity of 30 IU/ml (value obtained using **Magno Virus - Sacace** kit, [REF K-2-16/1000](#) starting from 1 ml of sample volume).

TROUBLESHOOTING


The results are not taken into account in the following cases:

1. If the Ct value obtained for the Positive Controls of Extraction or Amplification in the JOE channel is absent or is greater than the specified concentration range value, then it is necessary to repeat the test (from the DNA extraction stage for all samples in which HBV DNA was not found).
2. If the Ct value is obtained for the Negative Control of Extraction in the JOE channel and/or for the Negative Control of Amplification in the FAM and JOE channels, then it is necessary to repeat the test.
3. If the calculated concentrations of Positive Control-1 HBV and Positive Control-2 HBV exceed the range specified in the DataCard, then it is necessary to repeat the test (from the DNA extraction stage) for all samples.


EXPLANATION OF SYMBOLS


 REF Catalogue Number

 RUO For *Research Use Only*


 LOT Lot Number


 Expiration Date

 Contains reagents

 Caution!

 VER Version

 Manufacturer

 Temperature limitation

*SaCycler® is trademarks of Sacace Biotechnologies

*CFX-96™ and iQ5™ are trademarks of Bio-Rad Laboratories

*Rotor-Gene™ Technology is a registered trademark of Corbett Research

*MX3000P® and MX3005P® are trademarks of Stratagene

*Applied Biosystems® is trademarks of Applied Biosystems Corporation

* SmartCycler® is a registered trademark of Cepheid



Sacace Biotechnologies Srl

via Scalabrini, 44 – 22100 – Como – Italy Tel +390314892927 Fax +390314892926

mail: info@sacace.com web: www.sacace.com